

Derma Roller for Beard – A Healthier Beard Just One Product Away

Since the start of time, beards have been the pinnacle of manhood and elegance. Growing a beard does not only make a man look younger and cooler, but you can use it to cover up minor blemishes. In the fashion industry, derma rollers are efficient in improving skin but can you use them to stimulate a thick and healthy beard? The answer is yes!

Say Goodbye to Patchy Beard

[Derma roller for beard](#) works through a robust technology called micro-needling. Tiny harmless needles pierce through the skin to enhance blood flow to that area. Adequate blood circulation means a greater supply of essential nutrients. Thus, nurturing a healthier and stronger beard.

Skin-friendly Materials

This Derma roller is remarkable because it's designed from the best materials available in the industry, titanium and stainless steel. When you roll the stainless steel needles on the beard area, it boosts collagen production. Hair follicles take up this magic protein to stimulate growth and bless you with a full beard.

How to Use it?

If the bells of whistles of this product were intriguing enough, its ease of use would only enhance it.

Step 1 - Before using, disinfect the Derma Roller for Beard.

Step 2 - After washing your face, roll it on your neck area, right and left, 10-15 times.

Step 3 - Once done, clean and store your product in a safe place.

Precautions

This product must only be used in the areas where you wish to grow a beard. It is not recommended to use it on irritated skin or active acne, lips, or under the eyes

is not recommended. In case you accidentally drop your derma roller, be sure to replace it.