

Stimulate the Growth of a Beard with Derma Roller for Beard

What better way to step out of a boyish youth and enter manhood than growing a beard? Facial hair holds power to transform the shape of a face and amplify your charms. The presence of derma rollers as a tool to improve women's skin concerns is not unheard of. However, most have not yet been exposed to its magic in stimulating a thick and full beard in the men's industry.

[Derma rollers for Beard](#) function by using a robust technology called micro-needling. Covered in small titanium and stainless-steel needles, these rollers pierce through the skin's outermost layer. This procedure increases collagen, which is used up by hair follicles to promote the growth and quality of hair.

The skin underneath facial hair supports and dictates the quality of the beard. An adequate blood supply delivers active ingredients that strengthen the growth and quality of this hair. Derma rollers are known to be efficient in promoting blood circulation, and in return, this helps nurture a thick and manly beard.

The ease of use that this Derma Roller for beard provides makes this product so loveable. The length of the needles is a crucial factor in determining how frequently this product can be used. 0.5 mm needles are considered larger and should be limited to a maximum of two uses per week, whereas you should use smaller needles of 0.25mm more frequently. Roll 10-15 times on beard-growing areas and down your neck.

A few preventive measures can ensure the safety of this beard growth-stimulating tool. Wash your face and disinfect the roller before each use to prevent the spread of infection. Avoid rolling it on other parts like lips or under-eyes if you have active acne. If you accidentally drop this device, discard it and order a new one. It is also advised not to share it with anyone.